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Application

The touch screen of the CP9xx alarm indicator and operator panels enables the control of various operating theatre lights.

Some possible functions are: Switching surgical light(s) on/off, dimming, focusing, changing light intensity, changing light colour, camera control.

Operation is adapted to the functions of the respective operating theatre light. The scope of functions is specified by the respective manufacturer of the operating theatre light. Certain safety-critical functions can therefore only be performed on the light itself, such as controlling the laser function.

The scope of functions must be defined during project planning.

The specific visualisation elements of the individual manufacturers are integrated in the software of the alarm indicator and operator panel.

The list of manufacturers can be found in the table on the last page.

Particularities

Incorrect protocol messages between the alarm indicator and operator panel and the communication interface of the light manufacturers do not result in control of the light.

The lights can be controlled at any time via the control panel on the light.

Hardware

CP915, CP924 or a customised panel based on the CP9xx technology.

In some cases, a converter is required for communication with the manufacturer-specific interfaces of the operating theatre lights. Converters with galvanic separation must be used.

For the respective connection conditions, refer to the table on the last page with the approved and tested operating theatre lights.

In most cases the converter can be placed in the wall mounting enclosure of the CP9xx.

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The lights are connected in accordance with the respective manufacturer instructions.

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i *Caution! Currently, it is not possible to control operating theatre lights with a CP907 because the necessary USB port is missing.*

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i *Caution! In some cases, the converter cannot be placed in the panel enclosure, and this may result in the recommended cable length being exceeded.*

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Software

Communication is part of the COMTRAXX® software and is constantly maintained and developed. The respective surgical light type is implemented during project planning according to customer requirements.

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The integration of the operating theatre light and camera control into the visualisation of the alarm indicator and operator panel is carried out by partners of Bender GmbH & Co. KG.

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Adjustments are possible upon request.

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Visualisation

The user interface follows the design of the respective operating theatre light to make operation easier for the users.

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The respective user interface is part of the overall project planning and is placed after consultation with the customer.



Figure 1: Dräger Polaris 6xx

Figure 3: Merivaara Q-Flow

Figure 2: Trumple iLED7

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Components**Bender**

- Control panel, visualisation for control, RS232/USB adapter if required. Other adapters for interface conversion may be suitable

Operating theatre light manufacturers

- Operating theatre light(s), camera(s), communication interface, additional operating elements
- The operating theatre lights with all components are not maintained by Bender
- There are national exceptions

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Overview operating theatre lights

The list does not claim to be exhaustive. Throughout the product life cycle of alarm indicator and operator panels, new lights are constantly being integrated; others may no longer be available due to manufacturer discontinuations.

Hersteller Manufacturer	Name	Anzahl Leuchten Number OT lights	Kamera Camera	Schnittstelle Interface	Max. Länge [m] * Max. length [m] *
Dräger	Polaris 600	3	1	RS232 (G36180)	5
	Polaris Pro (+)	3	1	RS232 (GG10897)	5
KLS Martin	marLEDV	4	1	RS232	12
	marLED X	4	1	RS232	12
Getinge (Maquet)	Power LED 2	4	2	RS232	10
	Volista	4	2	RS232	10
	Volista Access II	4	2	RS232	10
Merivaara	Q-Flow	3	2	USB (xBee)	Wireless
Rimsa	Unica 520	3	1	RS232	12
	Unica 860	3	1	RS232	12
Simeon	Sim.LED 450 SC/MC	3	1	RS232	12
	Sim.LED 500 SC/MC	3	1	RS232	12
	Sim.LED 700 SC/MC	3	1	RS232	12
	Sim.LED 5000 SC/MC	3	1	RS232	12
	Sim.LED 7000 SC/MC	3	1	RS232	12
	Sim.LED 8000 SC/MC	3	1	RS232	12
Trumpf (Baxter Hillrom)	TruLight 5000	3	1	RS232	7.5
	iLED7	3	1	RS232	12

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All rights reserved.
Reprinting and duplicating only with
permission of the publisher.

© Bender GmbH & Co. KG, Germany
Subject to change! The specified
standards take into account the edition
valid until 07/2025 unless otherwise
indicated.