


# GRAND PROMERIX

Powerful and smooth Grand Promerix is heavy-duty, yet versatile operating table for all kind of surgical procedures. Improved ergonomic design improves the comfort of the whole surgical team. Comfort leads to concentration, confidence and better co-operation. That's efficiency.



## C-ARM ACCESS

The renewed slim base gives you perfect C-arm access. Grand Promerix meets even the strictest demands of modern operating theatres easily.



## IMPROVED ERGONOMICS

Grand Promerix has lots of improved details for better ergonomics, for example, hidden, smaller castors and inwards cutout for increased foot space. Human-centric design at its best.



### INCREASED HYGIENE

Details redesigned with patient safety in mind. Rethought frame make Grand Promerix hygienic and easy to clean.



### FLUENT LOCK™

Our Fluent locking system is both user-friendly and reliable. Connect with one click.



### HYGIENIC AND EASY TO CLEAN

New bellows and seat section make Grand Promerix hygienic and easy to clean.



### NEW DESIGN - IMPROVED DETAILS

Improved details improve hygiene and ergonomics for the operating room staff.

| Title                   | Sliding       |        | Non-Sliding   |        |
|-------------------------|---------------|--------|---------------|--------|
| Length (table top)      | 2100 mm       |        | 2100 mm       |        |
| Width (table top)       | 540 mm        |        | 540 mm        |        |
| Width (base)            | 652 mm        |        | 652 mm        |        |
| Height adjustment       | 650 - 1110 mm |        | 610 - 1070 mm |        |
| Max. patient weight     | 275 kg        | 600 lb | 325 kg        | 700 lb |
| Safe working load (SWL) | 275 kg        | 600 lb | 325 kg        | 700 lb |
| Trendelenburg           | -35 °         |        | -35 °         |        |
| Reverse Trendelenburg   | +35 °         |        | +35 °         |        |
| Head rest adjustment    | -45 - +45 °   |        | -45 - +45 °   |        |
| Back section adjustment | -45 - +80 °   |        | -45 - +80 °   |        |
| Leg section adjustment  | -105 - +70 °  |        | -95 - +70 °   |        |
| Castors size            | 80            |        | 80 mm         |        |
| Primary voltage         | 100 - 240 V   |        | 100 - 240 V   |        |
| Secondary voltage       | 24 V          |        | 24 V          |        |
| Max. power consumption  | 700 W         |        | 700 W         |        |
| Frequency               | 47-63 Hz      |        | 47-63 Hz      |        |
| Weight                  | 325 kg        |        | 300 kg        |        |
| Lateral tilt            | -25 - +25 °   |        | -20 - +20 °   |        |
| Lower frame length      | 1317 mm       |        | 1317 mm       |        |
| Max. lifting capacity   | 490 kg        |        | 490 kg        |        |
| Product Code            | 100060000     |        | 100060000     |        |